

# SENSOQUEST

**labcycler**

**Hightech Thermocycler**



**Purchase** SensoQuest Gradient Labcycler with Thermoblock 96 to receive **RM 1,500.00** worth of PCR Consumables\*

**Introductory Offer!**

Valid till 31 October 2011

## Winning Features:

- Attractive design, solid construction with 5.7" TFT touchscreen display
- Electroformed Gold plated silver thermoblocks
- Peltier elements tested to 600,000 cycles without failures
- Thermoblocks with 6-zone temperature regulation and self-calibrating temperature measuring system for precision thermal control
- Temperature range: -5°C to 99.9°C
- Automatic heated lid allows for fully programmable pressure and temperature, hotstart procedures compatible
- Intuitive user interface, software provides graphic monitoring of PCR process

## PCR Consumables from Axon Scientific

